




UVI HEALTH



Sirona

Shiatsu Neck and Back Massager
with Heating and Vibration

Important Safety Instructions

When using UVI Sirona make sure to follow all the instructions and necessary precautions listed in this user manual. UVI Sirona is unlikely to help with chronic pain, but it can have some therapeutic effects when used correctly. However, used in excess or inappropriately, it can cause injury. UVI Sirona should be used in moderation and should not be used for excessive periods of time.

Please read all safety instructions and warnings before using the product for the first time!



Do not operate UVI Sirona under a blanket or a pillow. Allow the heat to dissipate; otherwise the product can overheat, which can lead to injury, electric shock, or fire.



Do not use in environments with high temperatures or humidity since it could lead to electric shock or fire. If the product overheats, let it cool down for 30 minutes.




UVI Sirona should not be used by children or people with chronic illnesses, reduced physical, sensory or mental capabilities. Pregnant women should consult a doctor.





















When the product is not in use, please switch the power off. Do not leave the product unattended when plugged in, and never use the product while sleeping or napping.

UVI Sirona remote control



Product Operation Guide

1. To turn on the device, press the power button 
2. To adjust the height of the neck massage press  or .
3. To start the neck massage press . If you would like to change direction press  again. To turn off the neck massage press  for the third time.
4. Below you can find back massage functions.
Press once to turn on - to turn off press again:
 - To turn on the full back massage press 
 - To turn on the upper back massage press 
 - To turn on the lower back massage press 
 - To turn on the fixpoint kneading message press 
5. To turn on the regular massage mode press .
If you want to switch to waist massage press  again. If you press  for the third time, it will switch to back massage. To turn it off press  for the fourth time.
6. Press  to adjust timer to: 5 / 10 / 20 minutes.
7. To activate the vibration mode press . To switch between modes or turn it off press  again.
8. To turn on or off the heating function press .

Tips and Warnings

1. When using the product for the first time limit the use to 5 minutes or less, to allow your body to adjust to the experience.
2. The maximum recommended time to use the product is 30 minutes. After that, please wait for at least 15 minutes before using it again. Excessive use could lead to overheating or shortened lifespan of the product as well as injuries.
3. To soften the massage pressure, place a dry towel between yourself and the product.
4. Do not use the product if you have a swollen or inflamed area, a skin disorder, or an open wound.
5. UVI Sirona is not intended for medical use, nor shall it be used as a replacement for medical attention. The sole intended purpose of the product is to provide massage.

Cleaning and storage



Store UVI Sirona in a dry, dust-free location with with temperature between 15 do 30° and 30 - 70% relativ humidity. Ensure that the product is not exposed to direct sunlight.



Before cleaning, make sure that UVI Sirona has been turned off and that the adapter is plugged out from the outlet. Use a dry cloth to clean PU leather areas.



To prevent damage and discoloration, do not clean UVI Sirona with corrosive detergent. Do not wash it in the washing machine and keep it away from the water.



Do not apply excessive pressure to massage heads. If excessive force over 50 kilograms is applied directly, it may damage the product and invalidate the warranty.

If you encounter any problems, do not hesitate to contact us via our email info@uvihealth.com for help. More information and tips on how to use UVI Sirona to its fullest can be found at our website www.UVIHealth.com.

LEGIT d.o.o., Brnčičeva 13,
1231 Črnuče-Ljubljana, Slovenia

